

Employee of the month

The **H**igh **S**peed **S**orter –
My new high-speed colleague



Automatically more productive

Option to connect sample transport systems
(e.g. Tempus600®)

Bulk fill sample tubes in chutes
or the urgent sample
compartment – no inserting or
presorting

Gel-Check – prevents
recentrifugation and
minimises the risk of
incorrect analyses

Load and
unload
sample tubes
in the rack

Read barcode,
determine tube
type


Loop function –
automatic sample
feeding into the chute
of the Bulk Loader if no
order for the tubes has
been registered yet

HSS

BL

IO

ID

 SARSTEDT

→ SAMPLE LOADING

→ SAMPLE IDENTIFICATION/GEL-CHECK

The **High Speed Sorter** automates
your pre- and post-analytics

This creates time for other things.

Open screw caps and stoppers
regardless of manufacturer

Seal tubes for
transport and
long-term archiving

Sample sorting in racks
(closed or open)



DECAPPING

ALIQUOTING

RECAPPING

SORTING/DISTRIBUTING/ARCHIVING

Welcome to the team!



The automation solution for your workflow is just as personalised as your laboratory and your day-to-day requirements and challenges.

Every laboratory has to face new conditions and ask itself the following questions:

How do we want to deal with limited resources?

Which simple work processes can be automated and therefore optimised?


What sample throughput will this allow us to achieve?


Which functions are absolutely necessary for us?

Comparing this with our current situation: which analysis devices do we use and how much space do we have in our laboratory?

The modular design of the HSS (**H**igh **S**peed **S**orter) from SARSTEDT offers you the opportunity to create a solution precisely tailored to suit your needs.

Let's talk about your pre-analytics workflow.





- ✓ **Modular customisation** to suit the individual requirements of your laboratory environment
- ✓ **Space efficient** – optimal use of the available floorspace (can be positioned close to walls)
- ✓ **Tailor-made platform layouts** – tailored to your analysis devices

Automate your pre-analytics

- Sample loading
- Connection to a sample transport system (e.g. Tempus600®)
- Centrifugation status determination
- Sample identification
- Decapping
- Aliquoting
- Recapping
- Sorting
- Distributing
- Transfer to analysis platforms
- Archiving

Choose the optimal configuration for your laboratory

With the HSS system, SARSTEDT offers you customer-specific customisable automation solutions for your pre-analytics.

Take the strain off your employees

The HSS system allows relatively simple tasks and non-value creating activities, such as inserting and pre-sorting samples or opening and closing sample tubes, to be efficiently automated. This saves work for your employees who are (over)qualified for these tasks.

Optimise your processes


Heavy competition and cost pressure as well as limited personnel capacity means automating pre-analytical laboratory processes is particularly effective.

Learn more
about the
High Speed Sorter



hss.sarstedt.com

More than the sum of its parts



Bulk Loader module (BL)

Closed sample tubes are loaded manually into racks or automatically into chutes (Tempus600®) without pre-sorting or insertion.


The Bulk Loader has an integrated ID module for identifying samples and determining their centrifugation status (Gel-Check).

Automatic sample loading into the chute has also been integrated.



Decapper module (DC)

The Decapper removes screw caps and stoppers from various manufacturers – no pre-sorting, no cross-contamination.



Recapper module (RC)

The Recapper is available in 3 versions: it closes sample tubes with a diameter of 13 to 16 mm using a universal stopper, thus preventing evaporation. Alternatively, SARSTEDT tubes with a diameter of 13 to 15 mm are closed with the appropriate screw cap for transport or long-term archiving.



Plug & Play:
Automated pre- and
post-analytics with the
HSS modular system



Input/output module (IO)

The IO module allows open or closed sample tubes to be introduced into or ejected from the rack.

The platform layout is created specifically for each laboratory and, with its adapter plates (SARSTEDT FlexPlates), can be variably adjusted (e.g. to switch between routine and archive operation).



Identification module (ID)

The ID module identifies sample tubes from the barcode, determines the tube type (cap colour, shape and geometry) and establishes the centrifugation status (Gel-Check).



Sort Connect module (SC)

The Sort Connect module allows the HSS to be connected to an analysis platform (no need to manually transfer individual tubes or tubes in a rack).







Configure your HSS!

The modular concept allows for personalised configuration – tailor-made to suit your laboratory's needs.

Module options

Scope of functions, throughput, size, costs – the HSS laboratory automation system can also be individually resized retroactively.

Contact us about your individual configuration.





HSS – 6 modules

- Bulk packed sample loading
- Precise identification and sorting of non-centrifuged samples in centrifugation racks (Gel-Check)
- Loading after centrifugation in centrifugation or analysis racks
- Identification, selective opening and sorting
- Barcode alignment in analysis racks
- Selective closing with stoppers and archiving in racks
- Optional closing with screw caps



HSS – 5 modules

- Loading in racks for registration/recording or precise aliquoting downstream
- Identification, selective opening and sorting
- Sorting/distribution into rack types from various manufacturers
- Barcode alignment in analysis racks
- Selective closing with stoppers and archiving in racks
- Optional closing with screw caps



HSS – 4 modules

- Loading in various racks for registration/recording or precise aliquoting downstream
- Identification, selective opening and sorting
- Sorting/distribution into rack types from various manufacturers
- Barcode alignment in analysis racks
- Archiving of open samples possible



BL1200 – 4 modules

- Bulk packed sample loading
- Precise identification and pre-sorting in analysis racks (e.g. haematology, coagulation, urine)
- Barcode alignment in analysis racks
- Precise identification and sorting of gel samples in centrifugation racks with balancing tubes and buffer function
- Automated transfer of centrifugation racks or individual tubes to analysis platform systems



BL1200 – 3 modules

- Bulk packed sample loading
- Identification of centrifuged samples (Gel-Check)
- Selective opening and sorting in analysis racks
- Barcode alignment in analysis racks



BL1200 – 2 modules

- Bulk packed sample loading
- Sorting into centrifugation and analysis racks
- Barcode alignment in analysis racks
- Gel-Check optional

Technical details & data sheets



hss.sarstedt.com/technische-details

Even more efficient together: HSS and Tempus600®

Save time with the
Tempus600® sample
transport system



s-monovette-tempus600.
sarstedt.com



Tempus600® brings sample tubes to your laboratory simply, directly and securely.

Straight from the ward to the laboratory

With Tempus600® you accelerate sample transport, thus reducing the total turnaround time (ToTAT) for clinical samples and optimising the entire pre-analytical workflow.

This allows faster decisions to be made and can result in better treatment for patients.

Tempus600® reduces the gap from laboratory to ward.

Save with Tempus600®

- Save valuable time for your personnel with simple, direct sample transport and take the strain off existing transportation systems
- Use the potential of your HSS system even more efficiently
- Take the next step towards laboratory automation

Pre-analytical workflow made by SARSTEDT

SARSTEDT gives you holistic, coordinated products and systems that are consolidated in the pre-analytical workflow product range: from optimal sample quality, direct transport, automated sorting and distribution to reliable archiving. Utilise the synergies of our systems.

Let's work together.


Pre-analytical workflow made by SARSTEDT

Utilise the synergy of our
coordinated systems.

Discover the 360°
solutions for your
pre-analytics from
SARSTEDT



workflow.sarstedt.com



SARSTEDT AG & Co. KG

P.O. Box 12 20
D-51582 Nümbrecht

Phone: +49 2293 305 0
Fax: +49 2293 305 3992

Customer Service Germany
Phone +49 0800 0 83 305 0

export@sarstedt.com
www.sarstedt.com



More information,
videos and advice:
hss.sarstedt.com

If you have any questions:
We'd be happy to help!

Visit our website: www.sarstedt.com